

Notebook Explanation

1 Preparing the Dataset

Code

```
1 import glob
2 import argparse
3 import shutil
4
5
6 if __name__ == '__main__':
7     inp_folder = '/kaggle/input/alpha-dent/AlphaDent' + '/'
8     out_folder = '/kaggle/output/alpha-dent/AlphaDent_4_classes' + '/'
9
10    shutil.copytree(inp_folder, out_folder)
11
12    # Replace txt files
13    txt_paths = glob.glob(out_folder + '**/*.txt', recursive=True)
14    for txt_path in txt_paths:
15        lines = open(txt_path).readlines()
16        out = open(txt_path, 'w')
17        for line in lines:
18            if line[0] == '4' or line[0] == '5' or line[0] == '6' or
line[0] == '7' or line[0] == '8':
19                out.write('3' + line[1:])
20            else:
21                out.write(line)
22        out.close()
23
24    id_to_classes = {
25        1: 'Abrasion',
26        2: 'Filling',
27        3: 'Crown',
28        4: 'Caries',
29    }
30
31    # Create .yaml file
32    out = open(out_folder + 'yolo_seg_train.yaml', 'w')
33    out.write('path: {}\n'.format(out_folder))
34    out.write('train: images/train\n')
35    out.write('val: images/valid\n')
36    out.write('names:\n')
37    out.write('  0: Abrasion\n')
38    out.write('  1: Filling\n')
39    out.write('  2: Crown\n')
40    out.write('  3: Caries\n')
41    out.close()
42
43
44 print("Done")
```

Listing 1: Copying dataset and adjusting labels

Explanation

This section prepares the dataset for 4-class training by:

- Copying the original dataset into a new folder.
- Merging labels with IDs 4–8 into ID 3 (**Caries**).
- Writing a YOLO-compatible YAML file describing the dataset structure and class names.

2 Environment Setup and Dataset Splitting

Code and Explanation

1. Environment Setup

```
1 import os
2 import sys
3 import time
4 import glob
5 import cv2
6 import numpy as np
7 import pandas as pd
8 import matplotlib.pyplot as plt
9 import seaborn as sns
10 from pathlib import Path
11 from tqdm.auto import tqdm
12 import yaml
13 import random
14 from PIL import Image
15 import shutil
16 import warnings
17 warnings.filterwarnings('ignore')
18
19 # Disable wandb
20 os.environ['WANDB_DISABLED'] = 'true'
21
22 # Set random seeds for reproducibility
23 random.seed(42)
24 np.random.seed(42)
```

Listing 2: Import libraries, disable wandb, set random seeds

Explanation: This section imports all required libraries, disables `wandb`, and fixes random seeds to ensure reproducibility.

2. PyTorch and Ultralytics Setup

```
1 # Install required packages
2 print("Installing required packages...")
3 os.system('pip install -q ultralytics')
4
5 import torch
6 from ultralytics import YOLO
7
8 # Set deterministic behavior for PyTorch
9 torch.manual_seed(42)
10 if torch.cuda.is_available():
11     torch.cuda.manual_seed(42)
12     torch.backends.cudnn.deterministic = True
13     torch.backends.cudnn.benchmark = False
14
15 print(f'\nPyTorch Version: {torch.__version__}')
16 print(f'CUDA Available: {torch.cuda.is_available()}')
17 if torch.cuda.is_available():
18     print(f'CUDA Device: {torch.cuda.get_device_name(0)}')
```

Listing 3: Install ultralytics, set PyTorch deterministic behavior

Explanation: Installs Ultralytics YOLO, sets deterministic settings in PyTorch, and prints environment details.

3. Dataset Paths and Directory Creation

```
1 # Define original paths
2 BASE_PATH = '/kaggle/output/alpha-dent/AlphaDent_4_classes'
3 ORIGINAL_TRAIN_IMAGES_PATH = f'{BASE_PATH}/images/train'
4 VALID_IMAGES_PATH = f'{BASE_PATH}/images/valid'
5 TEST_IMAGES_PATH = f'{BASE_PATH}/images/test'
6 ORIGINAL_TRAIN_LABELS_PATH = f'{BASE_PATH}/labels/train'
7 VALID_LABELS_PATH = f'{BASE_PATH}/labels/valid'
8
9 # Output paths
10 OUTPUT_DIR = '/kaggle/working/'
11 WEIGHTS_DIR = f'{OUTPUT_DIR}/weights'
12 os.makedirs(WEIGHTS_DIR, exist_ok=True)
13
14 # Create new dataset structure with 90/10 split
15 NEW_DATASET_PATH = f'{OUTPUT_DIR}/alphadent_90_10_split'
16 NEW_TRAIN_IMAGES_PATH = f'{NEW_DATASET_PATH}/images/train'
17 NEW_EVAL_IMAGES_PATH = f'{NEW_DATASET_PATH}/images/eval'
18 NEW_VALID_IMAGES_PATH = f'{NEW_DATASET_PATH}/images/valid'
19 NEW_TEST_IMAGES_PATH = f'{NEW_DATASET_PATH}/images/test'
20 NEW_TRAIN_LABELS_PATH = f'{NEW_DATASET_PATH}/labels/train'
21 NEW_EVAL_LABELS_PATH = f'{NEW_DATASET_PATH}/labels/eval'
22 NEW_VALID_LABELS_PATH = f'{NEW_DATASET_PATH}/labels/valid'
23
24 # Create directories
25 os.makedirs(NEW_TRAIN_IMAGES_PATH, exist_ok=True)
26 os.makedirs(NEW_EVAL_IMAGES_PATH, exist_ok=True)
27 os.makedirs(NEW_VALID_IMAGES_PATH, exist_ok=True)
28 os.makedirs(NEW_TEST_IMAGES_PATH, exist_ok=True)
29 os.makedirs(NEW_TRAIN_LABELS_PATH, exist_ok=True)
30 os.makedirs(NEW_EVAL_LABELS_PATH, exist_ok=True)
31 os.makedirs(NEW_VALID_LABELS_PATH, exist_ok=True)
```

Listing 4: Define paths and create dataset structure

Explanation: Defines the original dataset structure, output directory, and creates new folders for the 90/10 split.

4. Train/Eval Split and File Copying

```
1 print("\n=== Creating 90/10 Train/Eval Split ===")
2
3 # Get all training images
4 original_train_images = sorted(glob.glob(f'{ORIGINAL_TRAIN_IMAGES_PATH}
5                               */*.jpg'))
6
7 print(f"Total original training images: {len(original_train_images)}")
8
9 # Shuffle and split into 90% train, 10% eval
10 random.shuffle(original_train_images)
11 split_idx = int(0.9 * len(original_train_images))
12 train_90_images = original_train_images[:split_idx]
13 eval_10_images = original_train_images[split_idx:]
14
15 print(f"90% for training: {len(train_90_images)}")
16 print(f"10% for evaluation: {len(eval_10_images)}")
17
18 def copy_files(image_list, dest_images_dir, dest_labels_dir,
19               source_labels_dir, desc):
20     """Copy images and corresponding labels to destination directories.
21     """
22     for img_path in tqdm(image_list, desc=desc):
23         img_filename = os.path.basename(img_path)
24         shutil.copy2(img_path, os.path.join(dest_images_dir,
25         img_filename))
26
27         label_filename = img_filename.replace('.jpg', '.txt')
28         source_label_path = os.path.join(source_labels_dir,
29         label_filename)
30         dest_label_path = os.path.join(dest_labels_dir, label_filename)
31
32         if os.path.exists(source_label_path):
33             shutil.copy2(source_label_path, dest_label_path)
```

Listing 5: Split training images into 90/10 and copy files

Explanation: Splits the training images into 90% train and 10% eval sets, then defines a helper function to copy images and their labels.

5. Copy Datasets

```
1 # Copy 90% training data
2 copy_files(train_90_images, NEW_TRAIN_IMAGES_PATH,
3            NEW_TRAIN_LABELS_PATH,
4            ORIGINAL_TRAIN_LABELS_PATH, "Copying 90% training data")
5
6 # Copy 10% evaluation data
7 copy_files(eval_10_images, NEW_EVAL_IMAGES_PATH, NEW_EVAL_LABELS_PATH,
8            ORIGINAL_TRAIN_LABELS_PATH, "Copying 10% evaluation data")
9
10 # Copy original validation data (unchanged)
11 original_valid_images = glob.glob(f'{VALID_IMAGES_PATH}/*.jpg')
12 copy_files(original_valid_images, NEW_VALID_IMAGES_PATH,
13            NEW_VALID_LABELS_PATH,
14            VALID_LABELS_PATH, "Copying validation data")
15
16 # Copy test data (images only, no labels)
17 print("Copying test images...")
18 test_images = glob.glob(f'{TEST_IMAGES_PATH}/*.jpg')
19 for img_path in tqdm(test_images, desc="Copying test images"):
20     img_filename = os.path.basename(img_path)
21     shutil.copy2(img_path, os.path.join(NEW_TEST_IMAGES_PATH,
22                                         img_filename))
```

Listing 6: Copy training, evaluation, validation, and test data

Explanation: Copies the datasets into the new structure: training, evaluation, validation, and test (images only).

6. YOLO Config Creation

```
1 # Define class information
2 CLASS_INFO = {
3     0: {'name': 'Abrasion', 'description': 'Teeth with mechanical wear
4         of hard tissues'},
5     1: {'name': 'Filling', 'description': 'Dental fillings of various
6         types'},
7     2: {'name': 'Crown', 'description': 'Dental crown (restoration)'},
8     3: {'name': 'Caries', 'description': 'Caries in fissures and pits'}
9 }
10
11 # Create YAML configuration for YOLO
12 print("\n=== Creating YOLO Configuration ===")
13 yolo_config = {
14     'path': NEW_DATASET_PATH,
15     'train': 'images/train',
16     'val': 'images/valid',
17     'test': 'images/test',
18     'nc': 4,
19     'names': [CLASS_INFO[i]['name'] for i in range(4)]
20 }
21
22 # Save the configuration
23 CUSTOM_YAML_PATH = f'{OUTPUT_DIR}/alphadent_config_90_10.yaml'
24 with open(CUSTOM_YAML_PATH, 'w') as f:
25     yaml.dump(yolo_config, f, default_flow_style=False)
26 print(f"Created custom YAML config at: {CUSTOM_YAML_PATH}")
```

Listing 7: Create YOLO config and dataset stats

Explanation: Defines 4 dental classes, creates a YOLO config file, and saves it in YAML format.

7. Class Distribution Analysis

```
1 def analyze_class_distribution(labels_path, dataset_name):
2     class_counts = {i: 0 for i in range(4)}
3     total_annotations = 0
4     label_files = glob.glob(f'{labels_path}/*.txt')
5
6     for label_file in tqdm(label_files, desc=f"Analyzing {dataset_name}
7     labels", leave=False):
8         if os.path.exists(label_file) and os.path.getsize(label_file) >
9         0:
10             with open(label_file, 'r') as f:
11                 for line in f:
12                     if line.strip():
13                         class_id = int(line.strip().split()[0])
14                         if 0 <= class_id < 4:
15                             class_counts[class_id] += 1
16                             total_annotations += 1
17
18     return class_counts, total_annotations
```

Listing 8: Analyze class distribution in datasets

Explanation: Defines a function to count class occurrences in label files, helping check dataset balance.

8. Training Configuration and Start

```
1 # Training configuration
2 print("\n=== Model Training Configuration ===")
3 EPOCHS = 30
4 IMAGE_SIZE = 640
5 BATCH_SIZE = 8 if torch.cuda.is_available() else 4
6 PATIENCE = 5
7
8 # Initialize and train model
9 print("\n=== Starting Model Training ===")
10 model = YOLO('yolov8x-seg.pt')
11
12 results = model.train(
13     data=CUSTOM_YAML_PATH,
14     epochs=EPOCHS,
15     imgsz=IMAGE_SIZE,
16     batch=BATCH_SIZE,
17     patience=PATIENCE,
18     save=True,
19     save_period=10,
20     project=OUTPUT_DIR,
21     name='alphadent_yolov8x_90_10',
22     exist_ok=True,
23     pretrained=True,
24     optimizer='AdamW',
25     lr0=0.001,
26     lrf=0.01,
27     momentum=0.937,
28     weight_decay=0.0005,
29     warmup_epochs=3.0,
30     warmup_momentum=0.8,
31     warmup_bias_lr=0.1,
32     box=7.5,
33     cls=0.5,
34     dfl=1.5,
35     hsv_h=0.015,
36     hsv_s=0.7,
37     hsv_v=0.4,
38     degrees=0.0,
39     translate=0.1,
40     scale=0.5,
41     shear=0.0,
42     perspective=0.0,
43     flipud=0.0,
44     fliplr=0.5,
45     mosaic=1.0,
46     mixup=0.0,
47     copy_paste=0.0,
48     plots=True,
49     device=0 if torch.cuda.is_available() else 'cpu',
50     workers=2,
51     verbose=True,
52     amp=True,
53     val=True
54 )
55
```

```
56 print("\nTraining completed!")
```

Listing 9: Training configuration and YOLOv8 training

Explanation: Defines hyperparameters (epochs, batch size, patience), initializes YOLOv8, and starts training with augmentation and optimization settings.

3 Model Loading, Validation and Evaluation

Code

```
1 # Load best model
2 print("=== Loading Best Model ===")
3 best_model_path = f'{OUTPUT_DIR}/alphadent_yolov8x_90_10/weights/best.
  pt'
4 if os.path.exists(best_model_path):
5     model = YOLO(best_model_path)
6     print(f"Loaded best model from: {best_model_path}")
7 else:
8     last_model_path = f'{OUTPUT_DIR}/alphadent_yolov8x_90_10/weights/
  last.pt'
9     if os.path.exists(last_model_path):
10        model = YOLO(last_model_path)
11        print(f"Loaded last model from: {last_model_path}")
12    else:
13        print("Warning: No trained model found, using pretrained model")
14    model = YOLO('yolov8x-seg.pt')
15
16 # Validate model on original validation set
17 print("\n=== Model Validation on Original Validation Set ===")
18 try:
19     metrics = model.val(
20         data=CUSTOM_YAML_PATH,
21         imgsz=IMAGE_SIZE,
22         batch=1,
23         conf=0.001,
24         iou=0.5,
25         max_det=300,
26         device=0 if torch.cuda.is_available() else 'cpu',
27         plots=False,
28         save_json=False,
29     )
30
31     print(f"\nValidation Results on Original Validation Set:")
32     print(f"mAP@50: {metrics.seg.map50:.4f}")
33     print(f"mAP@50-95: {metrics.seg.map:.4f}")
34 except Exception as e:
35     print(f"Validation error (non-critical): {e}")
36
37 # Create custom YAML for evaluation on the 10% held-out data
38 eval_config = {
39     'path': NEW_DATASET_PATH,
40     'train': 'images/train',
41     'val': 'images/eval', # Point to eval set for validation
42     'test': 'images/test',
43     'nc': 4,
```

```

44     'names': [CLASS_INFO[i]['name'] for i in range(4)]
45 }
46
47 EVAL_YAML_PATH = f'{OUTPUT_DIR}/alphadent_eval_config.yaml'
48 with open(EVAL_YAML_PATH, 'w') as f:
49     yaml.dump(eval_config, f, default_flow_style=False)
50
51 # Evaluate model on the 10% held-out labeled data
52 print("\n=== Model Evaluation on 10% Held-out Labeled Data ===")
53 try:
54     eval_metrics = model.val(
55         data=EVAL_YAML_PATH,
56         imsz=IMAGE_SIZE,
57         batch=1,
58         conf=0.001,
59         iou=0.5,
60         max_det=300,
61         device=0 if torch.cuda.is_available() else 'cpu',
62         plots=True,
63         save_json=True,
64         name='eval_10_percent'
65     )
66
67     print(f"\nEvaluation Results on 10% Held-out Data:")
68     print(f"mAP@50: {eval_metrics.seg.map50:.4f}")
69     print(f"mAP@50-95: {eval_metrics.seg.map:.4f}")
70     print(f"Per-class mAP@50:")
71     for i, map_val in enumerate(eval_metrics.seg.maps):
72         print(f"    {CLASS_INFO[i]['name']}: {map_val:.4f}")
73
74 except Exception as e:
75     print(f"Evaluation error: {e}")
76
77 # Final summary
78 print("\n=== Summary ===")
79 print(f"    Successfully created 90/10 train/eval split")
80 print(f"    Trained model on 90% of original training data ({len(
81     new_train_images)} images)")
81 print(f"    Evaluated model on 10% held-out labeled data ({len(
82     new_eval_images)} images)")
82 print(f"    Also validated on original validation set ({len(
83     new_valid_images)} images)")
83 print(f"    Model weights saved to: {OUTPUT_DIR}/
    alphadent_yolov8x_90_10/weights/")

```

Listing 10: Load best model, validate and evaluate on held-out data

Explanation

- Loads the best available checkpoint (best → last → pretrained) and runs validation.
- Validates on both the original validation set and the 10
- Writes a small eval YAML pointing to the held-out eval split.
- Prints a final summary of dataset sizes and where the model weights are stored.

4 Inference and Submission Pipeline

Code and Explanation

1. Clear Cache and Define Converter Function

```
1 import torch
2 torch.cuda.empty_cache()
3 # Inference on test set
4 print("\n=== Running Inference on Test Set ===")
5
6 def convert_to_submission_format(results, image_paths):
7     submission_rows = []
8
9     for idx, result in enumerate(results):
10         # Get image ID
11         image_id = os.path.basename(image_paths[idx]).replace('.jpg', '
12
13         if result.masks is not None and len(result.masks) > 0:
14             try:
15                 masks = result.masks.xy
16                 classes = result.bboxes.cls.cpu().numpy().astype(int)
17                 confidences = result.bboxes.conf.cpu().numpy()
18                 h, w = result.orig_shape
19
20                 # Process each detection
21                 for mask_idx in range(len(masks)):
22                     if mask_idx < len(classes) and mask_idx < len(
23 confidences):
24
25                         polygon = masks[mask_idx]
26
27                         if len(polygon) >= 3: # Valid polygon
28                             normalized_coords = []
29                             for point in polygon:
30                                 x_norm = float(point[0]) / w
31                                 y_norm = float(point[1]) / h
32                                 x_norm = max(0.0, min(1.0, x_norm))
33                                 y_norm = max(0.0, min(1.0, y_norm))
34                                 normalized_coords.extend([x_norm,
35 y_norm])
36
37                             poly_str = ' '.join([f'{coord:.6f}' for
38 coord in normalized_coords])
39
40                             submission_rows.append({
41                                 'patient_id': image_id,
42                                 'class_id': int(classes[mask_idx]),
43                                 'confidence': float(confidences[
44 mask_idx]),
45                                 'poly': poly_str
46                             })
47             except Exception as e:
48                 print(f"Error processing result for image {idx}: {e}")
49                 continue
```

```
46 return submission_rows
```

Listing 11: Empty CUDA cache and define conversion function

Explanation: First, GPU cache is cleared to avoid memory issues. Then, a function is defined to ****convert YOLO results into Kaggle submission format****: extracting polygons, normalizing coordinates to $[0, 1]$, and storing predictions as rows.

2. Run Inference on Test Images

```
1 # Process test images
2 test_images = sorted(glob.glob(f'{TEST_IMAGES_PATH}/*.jpg'))
3 all_submission_rows = []
4 INFERENCE_BATCH_SIZE = 2 if torch.cuda.is_available() else 4
5
6 print(f"Processing {len(test_images)} test images...")
7
8 for i in tqdm(range(0, len(test_images), INFERENCE_BATCH_SIZE)):
9     torch.cuda.empty_cache()
10    batch_images = test_images[i:i + INFERENCE_BATCH_SIZE]
11
12    try:
13        results = model.predict(
14            batch_images,
15            imgsz=IMAGE_SIZE,
16            conf=0.25,
17            iou=0.45,
18            max_det=300,
19            device=0 if torch.cuda.is_available() else 'cpu',
20            verbose=False,
21            agnostic_nms=True,
22            retina_masks=True,
23        )
24
25        batch_rows = convert_to_submission_format(results, batch_images)
26
27        all_submission_rows.extend(batch_rows)
28
29    except Exception as e:
30        print(f"Error in batch {i//INFERENCE_BATCH_SIZE}: {e}")
31        continue
32
33 print(f"\nGenerated {len(all_submission_rows)} predictions")
```

Listing 12: Process test set images in batches

Explanation: All test images are processed in ****small batches**** to avoid GPU memory overflow. The predictions are generated with YOLO, filtered by confidence (0.25) and IoU (0.45), then converted into submission rows using the function defined earlier.

3. Create Submission DataFrame

```
1 # Create submission DataFrame
2 print("\n=== Creating Submission File ===")
3 submission_df = pd.DataFrame(all_submission_rows)
4
5 # Ensure all test images are included
6 all_test_ids = [os.path.basename(img).replace('.jpg', '') for img in
7 test_images]
8 if len(submission_df) > 0:
9     predicted_ids = submission_df['patient_id'].unique()
10    missing_ids = set(all_test_ids) - set(predicted_ids)
11 else:
12    missing_ids = set(all_test_ids)
13
14 # Add dummy predictions if missing
15 if missing_ids:
16    print(f"Adding dummy predictions for {len(missing_ids)} images
17    without detections")
18    dummy_rows = []
19    for img_id in missing_ids:
20        dummy_rows.append({
21            'patient_id': img_id,
22            'class_id': 0,
23            'confidence': 0.01,
24            'poly': '0.1 0.1 0.1 0.2 0.2 0.2 0.2 0.1'
25        })
26
27    submission_df = pd.concat([submission_df, pd.DataFrame(dummy_rows)
28    ], ignore_index=True)
29
30 submission_df = submission_df.sort_values(['patient_id', 'confidence'],
31 ascending=[True, False])
32 submission_df = submission_df[['patient_id', 'class_id', 'confidence',
33 'poly']]
34 submission_df.to_csv('submission.csv', index=False)
35 print("Main submission file created: submission.csv")
```

Listing 13: Build submission file from predictions

Explanation: A DataFrame is created for submission. If some images have no detections, ****dummy polygons**** with low confidence are added to satisfy Kaggle format requirements. Finally, predictions are sorted and saved into **submission.csv**.

4. Verification of Submission

```
1 # Verify submission format
2 print("\n=== Verifying Submission Format ===")
3 print(f"Total predictions: {len(submission_df)}")
4 print(f"Unique images: {submission_df['patient_id'].nunique()}")
5 print(f"All test images included: {submission_df['patient_id'].nunique() == len(test_images)}")
6
7 print("\nFirst 5 rows of submission:")
8 print(submission_df.head())
9
10 # Check issues
11 print("\n=== Checking for Potential Issues ===")
12 missing_in_submission = set(all_test_ids) - set(submission_df['patient_id'].unique())
13 if missing_in_submission:
14     print(f"WARNING: Missing images in submission: {missing_in_submission}")
15 else:
16     print("    All test images have predictions")
17
18 # Class distribution
19 print("\nPredictions per class:")
20 class_dist = submission_df['class_id'].value_counts().sort_index()
21 for class_id, count in class_dist.items():
22     if 0 <= class_id < 9:
23         print(f"    Class {class_id} ({CLASS_INFO[class_id]['name']}): {count}")
24
25 # Confidence statistics
26 print(f"\nConfidence statistics:")
27 print(f"    Min: {submission_df['confidence'].min():.4f}")
28 print(f"    Max: {submission_df['confidence'].max():.4f}")
29 print(f"    Mean: {submission_df['confidence'].mean():.4f}")
30 print(f"    Median: {submission_df['confidence'].median():.4f}")
```

Listing 14: Verify submission consistency

Explanation: Ensures every test image has predictions, prints class distribution, and summarizes confidence values to check prediction quality.

5. High-Confidence Alternative Submission

```
1 # Create alternative submission with higher confidence threshold
2 print("\n=== Creating Alternative Submission (Higher Confidence) ===")
3 high_conf_df = submission_df[submission_df['confidence'] >= 0.3].copy()
4
5 # Ensure coverage of all images
6 high_conf_ids = high_conf_df['patient_id'].unique()
7 missing_high_conf = set(all_test_ids) - set(high_conf_ids)
8
9 if missing_high_conf:
10     for img_id in missing_high_conf:
11         img_preds = submission_df[submission_df['patient_id'] == img_id
12         ]
13         if len(img_preds) > 0:
14             high_conf_df = pd.concat([high_conf_df, img_preds.head(1)],
15             ignore_index=True)
16         else:
17             dummy_row = pd.DataFrame([
18                 'patient_id': img_id,
19                 'class_id': 0,
20                 'confidence': 0.01,
21                 'poly': '0.1 0.1 0.1 0.2 0.2 0.2 0.2 0.1'
22             ])
23             high_conf_df = pd.concat([high_conf_df, dummy_row],
24             ignore_index=True)
25
26 high_conf_df = high_conf_df.sort_values(['patient_id', 'confidence'],
27 ascending=[True, False])
28 high_conf_df.insert(0, 'id', range(1, len(high_conf_df) + 1))
29 high_conf_df.to_csv('submission_high_conf.csv', index=False)
30 print(f"Created high confidence submission with {len(high_conf_df)}
31 predictions")
32
33 print("\n=== Pipeline Completed Successfully! ===")
34 print("Submission files created:")
35 print(" - submission.csv (main submission)")
36 print(" - submission_high_conf.csv (alternative with higher confidence
37 threshold)")
```

Listing 15: Generate high-confidence version of submission

Explanation: An additional submission file is created with a ****confidence threshold of 0.3****. If some images are missing predictions at this threshold, the highest-confidence detection (or a dummy polygon) is added back. This produces `submission_high_conf.csv`, a stricter version of the submission.

5 Polygon-based Caries Classification Using ResNet18

This section describes the implementation of a polygon-based caries classifier built on top of a pre-trained ResNet18 model. The pipeline includes dataset preparation, model training, and evaluation.

5.1 Environment and Settings

We begin by importing the required libraries and setting paths, hyperparameters, and device configuration:

```
1 import os
2 import cv2
3 import torch
4 import numpy as np
5 from torch import nn
6 from torch.utils.data import Dataset, DataLoader
7 from torchvision import models, transforms
8 from PIL import Image
9 from glob import glob
10
11 # Paths and settings
12 TRAIN_PATH = "/kaggle/input/alpha-dent/AlphaDent/images/train"
13 TRAIN_LABEL_PATH = "/kaggle/input/alpha-dent/AlphaDent/labels/train"
14 NUM_CLASSES = 6 # caries1 -> caries6
15 IMG_SIZE = 224 # ResNet input size
16 BATCH_SIZE = 16
17 EPOCHS = 5
18 DEVICE = "cuda" if torch.cuda.is_available() else "cpu"
```

We define six caries sub-classes, set image size for ResNet input, batch size, and number of epochs. Computation runs on GPU if available.

5.2 Custom Dataset: Polygon Cropping

A custom Dataset class is implemented to load images and extract polygon regions corresponding to caries annotations:

```
1 class PolygonCariesDataset(Dataset):
2     def __init__(self, image_dir, label_dir, transform=None):
3         self.samples = []
4         self.transform = transform
5
6         image_paths = sorted(glob(os.path.join(image_dir, "*.jpg")))
7         label_paths = sorted(glob(os.path.join(label_dir, "*.txt")))
8
9         for img_path, lbl_path in zip(image_paths, label_paths):
10             with open(lbl_path, "r") as f:
11                 lines = f.readlines()
12                 for line in lines:
13                     parts = line.strip().split()
14                     class_id = int(parts[0])
15                     if 3 <= class_id <= 8: # merged caries classes
16                         coords = np.array(parts[1:], dtype=float).reshape
17                             (-1, 2)
18                         self.samples.append((img_path, coords, class_id-3))
19
20     def __len__(self):
21         return len(self.samples)
22
23     def __getitem__(self, idx):
24         img_path, coords, label = self.samples[idx]
25         img = cv2.imread(img_path)
26         img = cv2.cvtColor(img, cv2.COLOR_BGR2RGB)
```

```

26     h, w, _ = img.shape
27
28     # Crop polygon bounding box
29     xs, ys = coords[:,0], coords[:,1]
30     x1, y1 = max(int(xs.min()*w),0), max(int(ys.min()*h),0)
31     x2, y2 = min(int(xs.max()*w), w-1), min(int(ys.max()*h), h-1)
32     crop = img[y1:y2, x1:x2]
33
34     # Transform to tensor
35     crop = Image.fromarray(crop).convert("RGB")
36     if self.transform:
37         crop = self.transform(crop)
38     return crop, label

```

This class reads YOLO-format labels, filters only caries classes, crops bounding boxes from polygon masks, and returns cropped image patches with labels.

5.3 Data Transforms and Loader

The cropped patches are resized and normalized before being loaded into a `DataLoader`:

```

1 transform = transforms.Compose([
2     transforms.Resize((IMG_SIZE, IMG_SIZE)),
3     transforms.ToTensor(),
4     transforms.Normalize([0.5]*3, [0.5]*3)
5 ])
6
7 train_dataset = PolygonCariesDataset(TRAIN_PATH, TRAIN_LABEL_PATH,
8     transform=transform)
9 train_loader = DataLoader(train_dataset, batch_size=BATCH_SIZE, shuffle
10     =True)

```

5.4 Model: ResNet18 Fine-tuning

A pre-trained ResNet18 is adapted by replacing its final fully-connected layer with a classifier for six caries classes:

```

1 model = models.resnet18(pretrained=True)
2 model.fc = nn.Linear(model.fc.in_features, NUM_CLASSES)
3 model = model.to(DEVICE)

```

5.5 Training Setup

We use cross-entropy loss and Adam optimizer:

```

1 criterion = nn.CrossEntropyLoss()
2 optimizer = torch.optim.Adam(model.parameters(), lr=1e-4)

```

5.6 Training Loop

The model is trained for 5 epochs, reporting loss and accuracy per epoch:

```

1 for epoch in range(EPOCHS):
2     model.train()
3     running_loss, correct, total = 0.0, 0, 0
4

```

```

5     for imgs, labels in train_loader:
6         imgs, labels = imgs.to(DEVICE), labels.to(DEVICE)
7
8         optimizer.zero_grad()
9         outputs = model(imgs)
10        loss = criterion(outputs, labels)
11        loss.backward()
12        optimizer.step()
13
14        running_loss += loss.item()
15        _, predicted = outputs.max(1)
16        total += labels.size(0)
17        correct += (predicted == labels).sum().item()
18
19    print(f"Epoch [{epoch+1}/{EPOCHS}] "
20          f"Loss: {running_loss/len(train_loader):.4f} "
21          f"Accuracy: {correct/total:.4f}")

```

5.7 Saving the Model

After training, the model weights are saved:

```

1 torch.save(model.state_dict(), "caries_polygon_classifier.pth")

```

5.8 Validation and Metrics

The trained model is evaluated on a validation set, computing classification report and confusion matrix:

```

1 from sklearn.metrics import classification_report, confusion_matrix
2 import seaborn as sns
3 import matplotlib.pyplot as plt
4
5 val_dataset = PolygonCariesDataset(VAL_PATH, VAL_LABEL_PATH, transform=
    transform)
6 val_loader = DataLoader(val_dataset, batch_size=BATCH_SIZE, shuffle=
    False)
7
8 model.eval()
9 all_labels, all_preds = [], []
10 with torch.no_grad():
11     for imgs, labels in val_loader:
12         imgs, labels = imgs.to(DEVICE), labels.to(DEVICE)
13         outputs = model(imgs)
14         _, preds = outputs.max(1)
15         all_labels.extend(labels.cpu().numpy())
16         all_preds.extend(preds.cpu().numpy())
17
18 print(classification_report(all_labels, all_preds,
19                             target_names=[f"caries{i}" for i in range(1, NUM_CLASSES+1)]))
20
21 cm = confusion_matrix(all_labels, all_preds)
22 sns.heatmap(cm, annot=True, fmt="d", cmap="Blues",
23             xticklabels=[f"caries{i}" for i in range(1, NUM_CLASSES+1)
24             ],
25             yticklabels=[f"caries{i}" for i in range(1, NUM_CLASSES+1)
26             ])

```


```

25 plt.title("Confusion Matrix")
26 plt.show()

```


Table 1: Classification Report for Caries Classification

Class	Precision	Recall	F1-Score	Support
Caries1	0.74	0.82	0.78	62
Caries2	0.84	0.81	0.83	73
Caries3	0.78	0.64	0.70	33
Caries4	0.00	0.00	0.00	4
Caries5	0.78	0.86	0.82	81
Caries6	0.50	0.20	0.29	5
Accuracy			0.78	258
Macro Avg	0.61	0.56	0.57	258
Weighted Avg	0.77	0.78	0.77	258




6 Classification Examples


Prediction for test_043.jpg




Prediction for test_052.jpg



Prediction for test_040.jpg



Prediction for test_041.jpg



Prediction for test_014.jpg

